



# B 15

## Disinfection Wipes

### Product Information

#### Area of Application

B 15 Disinfection Wipes are alcohol-free and ready-to-use. They are particularly suitable for simultaneously disinfecting and cleaning surfaces of medical products (e.g. treatment units, medical equipment) in all fields in clinics, e.g. in operating rooms, intensive care units, casualty wards, central sterilization departments and for diagnosis devices which are also used for operations. B 15 Disinfection Wipes are also used by rescue services, in care homes (old people's homes, etc.), the food industry and in all fields where impeccable hygienic conditions are necessary.

#### Advantages of the product

- Alcohol-free ready-to-use wipes for simultaneously disinfecting and cleaning surfaces of medical products (e.g. treatment units, medical equipment, etc.) in clinics, practices and laboratories.
- Tested effectiveness and comprehensive effective spectrum against bacteria, fungi and all viruses (virucidal against enveloped and non-enveloped viruses).
- B 15 is VAH/DGHM-listed. Tested according to EN 13727, EN 13624, EN 14348, EN 14476, EN 14563.
- Extra large wipes: 175 x 280 mm.
- Intensive cleaning power and excellent material compatibility, even in case of alcohol-sensitive surfaces, e.g. acrylic glass.
- Practical and easy application, removal from dispenser box.
- Increased safety because of remnant effect.
- For easy disinfection of surfaces.

#### Microbiological activity

B 15 have the following characteristics of performance:

Application	Time
Surface disinfection (in accordance with VAH/DGHM)	2 min.
Bacteria incl. MRSA <sup>1)</sup> , fungi <sup>2)</sup>	2 min.
Tb bacteria	60 min.
Vaccinia viruses inclusive of HBV, HCV and HIV <sup>3)</sup>	1 min.
Adeno viruses <sup>4)</sup>	15 min.
Polyoma SV 40 viruses	5 min.
Polio viruses <sup>4)</sup>	5 min.
Noro viruses <sup>5)</sup>	5 min.

<sup>1)</sup> Effectiveness against bacteria implies effectiveness against MRSA.

<sup>2)</sup> Test for *Candida albicans* and *Aspergillus niger*.

<sup>3)</sup> In accordance with RKI recommendations (Bundesgesundheitsbl. 47, 62- 66, 2004).

<sup>4)</sup> Values without organic load.

<sup>5)</sup> Effectiveness against poliovirus implies effectiveness against norovirus.

#### Application

Moisten surface and allow product to take effect for 2 minutes in accordance with VAH/DGHM. We recommend 15 minutes for inactivating all viruses. Evenly moisten surfaces to be treated.

#### Notes

Do not use B 15 Disinfection Wipes on surfaces made of aluminium.

#### Chemical composition

B 15 Disinfection Wipes contain B 15. 1 wipe contains 0.029 g 3-aminopropyl-dodecyl-1,3-propane diamine, 0.017 g didecyl-methyl-poly(oxyethyl)-ammonium propionate, non-ionic surfactants, alkaline cleaning components, complexing agents, adjuvants and water.

#### Packing units

Basic set: reusable dispenser box containing 1 x 50 wipes  
refill pack: carton of 4 x 50 wipes

#### Shelf life

4 years

#### Please turn over

Storage • Physical properties • Toxicology • Ecology • Caution



### Storage

Store the preparation in a cool place, but not below 5°C.

### Physical properties

#### Solution of active ingredients:

Appearance: clear solution  
pH: 11.6 ± 0.3

### Toxicology

Based on the data on hand for B 15, the B 15 Disinfection Wipes are classified as "non-toxic". There is no risk from inhalation under normal conditions of use.

### Ecology

B 15 Disinfection Wipes are ecologically safe, as the organic ingredients contained in B 15 are biodegradable when appropriately diluted in water. The packing is made of polyethylene and polyester and can easily be exploited thermally. Dispose of used wipes in normal domestic waster.

### Caution

Avoid contact with skin and eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. After contact with skin, wash immediately with plenty of water. Wear suitable gloves.

### CE 0124

The information on our products is based on detailed internal and external scientific examinations, and recommendations are given to the best of our knowledge. Liability is being judged in accordance with our General Terms of Sale and Delivery which are available at any time upon request.

### Complementary products from oro<sup>®</sup> Hygienesystem

Instrument disinfection • Disinfection of surfaces and inventory • Disinfection and care of hands and skin • Disinfection of special applications • Accessories

Information as of 01/10